

CD63 Rabbit mAb

ABZ2878


Product Overview

Name	CD63 Rabbit mAb
Accession(Primary)	P08962
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 26 kDa; Observed MW: 26-65 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:50-1:100


Target information (P08962)

Gene Name CD63

About Gene Symbol



Western blot analysis of CD63 in HeLa lysates using CD63 antibody.



Immunohistochemistry analysis of paraffin-embedded Human breast cancer using CD63 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.