

CD7 Rabbit mAb


ABZ2879

Product Overview


Name	CD7 Rabbit mAb
Accession(Primary)	P09564
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 25 kDa; Observed MW: 35-40 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:1000-1:5000; ICC/IF-1:20-1:50; FC-1:50-1:100; IP-1:10-1:100

Target information (P09564)


Gene Name	CD7
About Gene Symbol	



Immunocytochemistry analysis of CD7 (green) in Jurkat using CD7 antibody, and DAPI (blue).



Western blot analysis of CD7 in Jurkat lysates using CD7 antibody.



Immunohistochemistry analysis of paraffin-embedded Human tonsil using CD7 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.