


CDC7 Rabbit mAb

ABZ2883

Product Overview

Name	CDC7 Rabbit mAb
Accession(Primary)	O00311
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 64 kDa; Observed MW: 64 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:200-1:500; IP-1:20-1:50

Target information (O00311)

Gene Name CDC7**About Gene Symbol**

Western blot analysis of CDC7 in HeLa, CHO-K1, Jurkat lysates using CDC7 antibody.

