

CDK1 Rabbit mAb


ABZ2884

Product Overview

Name	CDK1 Rabbit mAb
Accession(Primary)	P06493
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 34 kDa; Observed MW: 34 kDa
Recommended Dilution	WB-1:1000-1:50000; IHC-1:50-1:2000; ICC/IF-1:100-1:200; FC-1:50-1:100; IP-1:10-1:100

Target information (P06493)

Gene Name	CDK1
About Gene Symbol	



Western blot analysis of CDK1 in K562, HeLa
lysates using CDK1 antibody.