

CDK4 Rabbit mAb


ABZ2886

Product Overview


Name	CDK4 Rabbit mAb
Accession(Primary)	P11802
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 34 kDa; Observed MW: 34 kDa
Recommended Dilution	WB-1:500-1:1000; ICC/IF-1:100-1:200; IP-1:20-1:50; IHC-1:50-1:100

Target information (P11802)

Gene Name	CDK4
About Gene Symbol	CDK4 (Cyclin Dependent Kinase 4) is a protein kinase and signaling regulator involved in cell cycle regulation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of Cdk4 in K562, C6, 3T3, HeLa lysates using CDK4 antibody.



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using Cdk4 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.