

Chk1 Rabbit mAb


ABZ2891

Product Overview


Name	Chk1 Rabbit mAb
Accession(Primary)	O14757
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 54 kDa; Observed MW: 54 kDa
Recommended Dilution	WB-1:500-1:1000; ICC/IF-1:20-1:50; FC-1:20; IP-1:20

Target information (O14757)

Gene Name	CHEK1
About Gene Symbol	CHEK1 (Checkpoint Kinase 1) is a serine/threonine protein kinase involved in cell cycle regulation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of Chk1 in K562 lysates using Chk1 antibody.



Western blot analysis of Chk1 in K562, 3T3, HeLa lysates using Chk1 antibody.