

Chromogranin A Rabbit mAb


ABZ2893

Product Overview

Name	Chromogranin A Rabbit mAb
Accession(Primary)	P10645
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 51 kDa; Observed MW: 92 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:100-1:200; ICC/IF-1:100-1:200; FC-1:100-1:200

Target information (P10645)

Gene Name	CHGA
About Gene Symbol	CHGA (Chromogranin A) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of Chromogranin A in PC-12
lysates using Chromogranin A antibody.