

## Claudin 3 Rabbit mAb


ABZ2897

### Product Overview

|                             |                                                                                      |
|-----------------------------|--------------------------------------------------------------------------------------|
| <b>Name</b>                 | Claudin 3 Rabbit mAb                                                                 |
| <b>Accession(Primary)</b>   | O15551                                                                               |
| <b>Host Species</b>         | Rabbit                                                                               |
| <b>Clonality</b>            | Monoclonal                                                                           |
| <b>Isotype</b>              | IgG                                                                                  |
| <b>Species Reactivity</b>   | Human                                                                                |
| <b>Calculated MW</b>        | Calculated MW: 23 kDa; Observed MW: 20 kDa                                           |
| <b>Recommended Dilution</b> | WB-1:500-1:1000; IHC-P-1:100-1:200;<br>ICC/IF-1:50-1:100; FC-1:20-1:50; IP-1:20-1:50 |

### Target information (O15551)

|                          |       |
|--------------------------|-------|
| <b>Gene Name</b>         | CLDN3 |
| <b>About Gene Symbol</b> |       |



Western blot analysis of Claudin 3 in MCF-7  
lysates using Claudin 3 antibody.