

Clusterin alpha Chain Rabbit mAb

ABZ2899

Product Overview

Name	Clusterin alpha Chain Rabbit mAb
Accession(Primary)	Q06890
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Mouse, Rat
Calculated MW	Calculated MW: 52 kDa; Observed MW: 42 kDa
Recommended Dilution	WB-1:500-1:1000; IP-1:20; IHC-1:50-1:100

Target information (Q06890)

Gene Name

Clu

About Gene Symbol



