


COG6 Rabbit mAb

ABZ2900

Product Overview

Name	COG6 Rabbit mAb
Accession(Primary)	Q9Y2V7
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 73 kDa; Observed MW: 73 kDa
Recommended Dilution	WB-1:1000-1:5000

Target information (Q9Y2V7)

Gene Name COG6**About Gene Symbol**

Western blot analysis of COG6 in K562, 3T3
lysates using COG6 antibody.

