

CREB Rabbit mAb


ABZ2903

Product Overview


Name	CREB Rabbit mAb
Accession(Primary)	P16220
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 35 kDa; Observed MW: 43 kDa
Recommended Dilution	WB-1:500-1:10000; IHC-1:50-1:2000; ICC/IF-1:100-1:200; FC-1:50-1:100; IP-1:20-1:100

Target information (P16220)


Gene Name	CREB1
About Gene Symbol	CREB1 (CAMP Responsive Element Binding Protein 1) is a transcription factor involved in gene transcription and regulation. It is an important research target in gene regulation, development, and cancer research and is widely studied to elucidate cellular



Western blot analysis of CREB in 3T3 lysates using CREB antibody.



Western blot analysis of CREB in HeLa lysates using CREB antibody.



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using CREB antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.