

CYB5R3 Rabbit mAb

ABZ2908


Product Overview


Name	CYB5R3 Rabbit mAb
Accession(Primary)	P00387
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Rat
Calculated MW	Calculated MW: 34 kDa; Observed MW: 34 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:200-1:500

Target information (P00387)

Gene Name CYB5R3

About Gene Symbol





Immunohistochemistry analysis of paraffin-embedded Human breast cancer using CYB5R3 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.