

Cyclin A2 Rabbit mAb


ABZ2909

Product Overview


Name	Cyclin A2 Rabbit mAb
Accession(Primary)	P20248
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 49 kDa; Observed MW: 49 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:200-1:500; ICC/IF-1:100-1:200; IP-1:20-1:50

Target information (P20248)


Gene Name	CCNA2
About Gene Symbol	CCNA2 (Cyclin A2) is a protein kinase and signaling regulator involved in cell cycle regulation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of Cyclin A2 in K562, C6, 3T3, HeLa lysates using Cyclin A2 antibody.



Immunocytochemistry analysis of Cyclin A2(green) in HeLa using Cyclin A2 antibody, and DAPI(blue)



Immunohistochemistry analysis of paraffin-embedded Human tonsil using Cyclin A2 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.