

Cyclophilin B Rabbit mAb

ABZ2911


Product Overview

Name	Cyclophilin B Rabbit mAb
Accession(Primary)	P23284
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 24 kDa; Observed MW: 20 kDa
Recommended Dilution	WB-1:1000-1:5000


Target information (P23284)

Gene Name	PPIB
About Gene Symbol	PPIB (Peptidylprolyl Isomerase B) is a peptidylprolyl isomerase (cyclophilin) involved in cellular metabolism. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.

Mouse pancreas



Western blot analysis of Cyclophilin B in mouse pancreas lysates using Cyclophilin B antibody.



Western blot analysis of Cyclophilin B in HeLa, C6, rat Brain lysates using Cyclophilin B antibody.