

Cytochrome P450 17A1/CYP17A1 Rabbit mAb


ABZ2913

Product Overview


Name	Cytochrome P450 17A1/CYP17A1 Rabbit mAb
Accession(Primary)	P05093
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 57 kDa; Observed MW: 57 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC:1:20; ICC/IF-1:20; FC-1:20; IP-1:20

Target information (P05093)


Gene Name	CYP17A1
About Gene Symbol	



Immunohistochemistry analysis of paraffin-embedded mouse kidney using Cytochrome P450 17A1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Western blot analysis of Cytochrome P450 17A1/CYP17A1 in K562, C6, 3T3 lysates using Cytochrome P450 17A1/CYP17A1 antibody.



Immunocytochemistry analysis of Cytochrome P450 17A1(green) in HeLa cells using Cytochrome P450 17A1 antibody, and DAPI(blue)