

Cytokeratin 18 Rabbit mAb


ABZ2914

Product Overview


Name	Cytokeratin 18 Rabbit mAb
Accession(Primary)	P05783
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 48 kDa; Observed MW: 48 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-1:100-1:500; IF-1:50-1:200

Target information (P05783)


Gene Name	KRT18
About Gene Symbol	



Western blot analysis of Cytokeratin 18 in HeLa lysates using Cytokeratin 18 antibody.



Immunohistochemistry analysis of paraffin-embedded Human Cholangiocarcinoma using Cytokeratin 18 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunocytochemistry analysis of Cytokeratin 18 (green) in hela using Cytokeratin 18 antibody, and DAPI(blue)