

Cytokeratin 19 Rabbit mAb

ABZ2915


Product Overview

Name	Cytokeratin 19 Rabbit mAb
Accession(Primary)	P08727
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 44 kDa; Observed MW: 41 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC:1:20; ICC/IF-1:20


Target information (P08727)

Gene Name KRT19

About Gene Symbol



Western blot analysis of Cytokeratin 19 in K562 lysates using Cytokeratin 19 antibody.



Immunohistochemistry analysis of paraffin-embedded Human colon cancer using Cytokeratin 19 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.