

# Cytokeratin 19 Rabbit mAb


ABZ2916

## Product Overview


|                             |                                                                                    |
|-----------------------------|------------------------------------------------------------------------------------|
| <b>Name</b>                 | Cytokeratin 19 Rabbit mAb                                                          |
| <b>Accession(Primary)</b>   | P08727                                                                             |
| <b>Host Species</b>         | Rabbit                                                                             |
| <b>Clonality</b>            | Monoclonal                                                                         |
| <b>Isotype</b>              | IgG                                                                                |
| <b>Species Reactivity</b>   | Human                                                                              |
| <b>Calculated MW</b>        | Calculated MW: 44 kDa; Observed MW: 41 kDa                                         |
| <b>Recommended Dilution</b> | WB-1:1000-1:5000; IHC-P-1:200-1:1000;<br>ICC/IF-1:100-1:200; FC-1:20-1:50; IP-1:20 |

## Target information (P08727)


|                          |       |
|--------------------------|-------|
| <b>Gene Name</b>         | KRT19 |
| <b>About Gene Symbol</b> |       |



Western blot analysis of Cytokeratin 19 in rat Brain lysates using Cytokeratin 19 antibody.



Immunohistochemistry analysis of paraffin-embedded Human colon cancer using Cytokeratin 19 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunocytochemistry analysis of Cytokeratin 19 (green) in MCF-7 using Cytokeratin 19 antibody, and DAPI(blue)