

DDAH1 Rabbit mAb


ABZ2919

Product Overview

Name	DDAH1 Rabbit mAb
Accession(Primary)	O94760
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 31 kDa; Observed MW: 31 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:50; ICC/IF-1:50-1:100; FC-1:20-1:50; IP-1:20-1:50

Target information (O94760)

Gene Name	DDAH1
About Gene Symbol	DDAH1 (Dimethylarginine Dimethylaminohydrolase 1) is a metabolic enzyme involved in cellular metabolism. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of DDAH1 in HeLa lysates
using DDAH1 antibody.