

## delta 1 Catenin Rabbit mAb


ABZ2920

### Product Overview


|                             |                                                                   |
|-----------------------------|-------------------------------------------------------------------|
| <b>Name</b>                 | delta 1 Catenin Rabbit mAb                                        |
| <b>Accession(Primary)</b>   | O60716                                                            |
| <b>Host Species</b>         | Rabbit                                                            |
| <b>Clonality</b>            | Monoclonal                                                        |
| <b>Isotype</b>              | IgG                                                               |
| <b>Species Reactivity</b>   | Human, Mouse                                                      |
| <b>Calculated MW</b>        | Calculated MW: 108 kDa; Observed MW: 108 kDa                      |
| <b>Recommended Dilution</b> | WB-1:1000-1:5000; IHC-P-1:20-1:50;<br>ICC/IF-1:100-1:200; FC-1:20 |

### Target information (O60716)

|                          |                                                                                                                                                                                                                                                                                               |
|--------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Gene Name</b>         | CTNND1                                                                                                                                                                                                                                                                                        |
| <b>About Gene Symbol</b> | CTNND1 (Catenin Delta 1) is a cell junction/adhesion protein involved in cytoskeletal organization. It is an important research target in cell biology, cancer metastasis, and cytoskeletal research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets. |



Western blot analysis of delta 1 Catenin in C6, 3T3, HeLa lysates using delta 1 Catenin antibody.



Immunohistochemistry analysis of paraffin-embedded Human colon cancer using delta 1 Catenin antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.