

DFNA5/GSDME Rabbit mAb


ABZ2922

Product Overview

Name	DFNA5/GSDME Rabbit mAb
Accession(Primary)	O60443
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 55 kDa; Observed MW: 55 kDa
Recommended Dilution	WB-1:500-1:1000; FC-1:50-1:100; IP-1:20-1:50

Target information (O60443)

Gene Name	GSDME
About Gene Symbol	GSDME (Gasdermin E) is a tumor suppressor protein involved in receptor-mediated signaling. It is an important research target in cancer research and hearing loss and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of DFNA5/GSDME in U251MG, U2OS, HepG2, HeLa lysates using GSDME antibody.