

DHRS1 Rabbit mAb

ABZ2924


Product Overview

Name	DHRS1 Rabbit mAb
Accession(Primary)	Q96LJ7
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 34 kDa; Observed MW: 34 kDa
Recommended Dilution	WB-1:500-1:1000

Target information (Q96LJ7)

Gene Name DHRS1

About Gene Symbol



Western blot analysis of DHRS1 in C6 lysates using DHRS1 antibody.

