


DHX8 Rabbit mAb

ABZ2925

Product Overview

Name	DHX8 Rabbit mAb
Accession(Primary)	Q14562
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 139 kDa; Observed MW: 139 kDa
Recommended Dilution	WB-1:500-1:1000; IF-1:50-1:200

Target information (Q14562)

Gene Name DHX8**About Gene Symbol**

Western blot analysis of DHX8 in HeLa, 293T lysates using DHX8 antibody.

