


DjC8 Rabbit mAb

ABZ2926


Product Overview

Name	DjC8 Rabbit mAb
Accession(Primary)	O75937
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 30 kDa; Observed MW: 30 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:50-1:100

Target information (O75937)

Gene Name DNAJC8**About Gene Symbol**

Western blot analysis of DjC8 in mouse brain, rat lung lysates using DNAJC8 antibody.



Western blot analysis of DjC8 in Jurkat, HeLa, COS7, Raji lysates using DjC8 antibody.