


DNAJA4 Rabbit mAb


ABZ2927

Product Overview

Name	DNAJA4 Rabbit mAb
Accession(Primary)	Q8WW22
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 45 kDa; Observed MW: 45 kDa
Recommended Dilution	WB-1:500-1:1000; IP-1:20

Target information (Q8WW22)

Gene Name DNAJA4**About Gene Symbol**



Western blot analysis of DNAJA4 in C6, rat Brain lysates using DNAJA4 antibody.