

DNAJB6 Rabbit mAb


ABZ2928

Product Overview


Name	DNAJB6 Rabbit mAb
Accession(Primary)	O75190
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 36 kDa; Observed MW: 36 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:50-1:100; ICC/IF-1:100-1:200; FC-1:200-1:500; IP-1:100-1:200

Target information (O75190)

Gene Name	DNAJB6
About Gene Symbol	



Immunohistochemistry analysis of paraffin-embedded mouse colon using DNAJB6 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Western blot analysis of DNAJB6 in K562, rat Brain, C6, 3T3, HeLa lysates using DNAJB6 antibody.