


DOCK2 Rabbit mAb

ABZ2929

Product Overview

Name	DOCK2 Rabbit mAb
Accession(Primary)	Q92608
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 212 kDa; Observed MW: 212 kDa
Recommended Dilution	WB-1:1000-1:5000

Target information (Q92608)

Gene Name DOCK2**About Gene Symbol**

Western blot analysis of DOCK2 in Jurkat lysates using DOCK2 antibody.

