

DOCK8 Rabbit mAb


ABZ2930

Product Overview

Name	DOCK8 Rabbit mAb
Accession(Primary)	Q8NF50
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 239 kDa; Observed MW: 239 kDa
Recommended Dilution	WB-1:1000-1:5000; ICC/IF-1:200-1:500; FC-1:10-1:100

Target information (Q8NF50)

Gene Name	DOCK8
About Gene Symbol	



Western blot analysis of DOCK8 in Ramos lysates using DOCK8 antibody.