

DOHH Rabbit mAb

ABZ2931


Product Overview

Name	DOHH Rabbit mAb
Accession(Primary)	Q9BU89
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 33 kDa; Observed MW: 33 kDa
Recommended Dilution	WB-1:500-1:1000; IP-1:20

Target information (Q9BU89)

Gene Name DOHH

About Gene Symbol



Western blot analysis of DOHH in K562, HeLa lysates using DOHH antibody.

