

eIF4E Rabbit mAb


ABZ2932

Product Overview


Name	eIF4E Rabbit mAb
Accession(Primary)	P06730
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 25 kDa; Observed MW: 25 kDa
Recommended Dilution	WB-1:500-1:1000; IHC:1:50-1:100; ICC/IF-1:200-1:500; FC-1:100-1:200; IP-1:20

Target information (P06730)


Gene Name	EIF4E
About Gene Symbol	EIF4E (Eukaryotic Translation Initiation Factor 4E) is a protein kinase and signaling regulator involved in translational regulation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Immunocytochemistry analysis of eIF4E (green) in HCT116 using eIF4E antibody, and DAPI(blue).



Western blot analysis of eIF4E in HL-60 lysates using eIF4E antibody.



Western blot analysis of eIF4E in C6, 3T3 lysates
using eIF4E antibody