

## eIF4E Rabbit mAb


ABZ2933

### Product Overview


|                             |                                                                              |
|-----------------------------|------------------------------------------------------------------------------|
| <b>Name</b>                 | eIF4E Rabbit mAb                                                             |
| <b>Accession(Primary)</b>   | P06730                                                                       |
| <b>Host Species</b>         | Rabbit                                                                       |
| <b>Clonality</b>            | Monoclonal                                                                   |
| <b>Isotype</b>              | IgG                                                                          |
| <b>Species Reactivity</b>   | Human, Mouse, Rat                                                            |
| <b>Calculated MW</b>        | Calculated MW: 25 kDa; Observed MW: 25 kDa                                   |
| <b>Recommended Dilution</b> | WB-1:500-1:1000; IHC:1:200-1:500;<br>ICC/IF-1:20-1:50; FC-1:20-1:50; IP-1:20 |

### Target information (P06730)

|                          |                                                                                                                                                                                                                                                                                                                              |
|--------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Gene Name</b>         | EIF4E                                                                                                                                                                                                                                                                                                                        |
| <b>About Gene Symbol</b> | EIF4E (Eukaryotic Translation Initiation Factor 4E) is a protein kinase and signaling regulator involved in translational regulation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets. |



Immunohistochemistry analysis of paraffin-embedded mouse colon using eIF4E antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Western blot analysis of eIF4E in K562, rat Brain, C6, 3T3, HeLa lysates using eIF4E antibody.