

ELK1 Rabbit mAb


ABZ2934

Product Overview


Name	ELK1 Rabbit mAb
Accession(Primary)	P19419
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 45 kDa; Observed MW: 62 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:20-1:50; ICC/IF-1:200-1:500; FC-1:20-1:50

Target information (P19419)

Gene Name	ELK1
About Gene Symbol	ELK1 (ETS Transcription Factor ELK1) is a receptor tyrosine kinase involved in cell signaling and proliferation. It is an important research target in cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of ELK1 in HeLa lysates using ELK1 antibody.



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using ELK1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.