

ELK1 Rabbit mAb


ABZ2935

Product Overview

Name	ELK1 Rabbit mAb
Accession(Primary)	P19419
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 45 kDa; Observed MW: 62 kDa
Recommended Dilution	WB-1:500-1:1000; ICC/IF-1:50-1:100

Target information (P19419)

Gene Name	ELK1
About Gene Symbol	ELK1 (ETS Transcription Factor ELK1) is a receptor tyrosine kinase involved in cell signaling and proliferation. It is an important research target in cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of ELK1 in C6, 3T3, HeLa
lysates using ELK1 antibody.