

EN1 Rabbit mAb

ABZ2936


Product Overview

Name	EN1 Rabbit mAb
Accession(Primary)	Q05925
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 40 kDa; Observed MW: 50 kDa
Recommended Dilution	WB-1:1000-1:5000

Target information (Q05925)

Gene Name EN1

About Gene Symbol



Western blot analysis of EN1/Engrailed 1 in rat Brain lysates using EN1 antibody.

