

ENO3 Rabbit mAb


ABZ2937

Product Overview


Name	ENO3 Rabbit mAb
Accession(Primary)	P13929
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 47 kDa; Observed MW: 47 kDa
Recommended Dilution	WB-1:500-1:1000; ICC/IF-1:200-1:500; FC-1:50-1:100

Target information (P13929)

Gene Name	ENO3
About Gene Symbol	



Western blot analysis of ENO3 in rat Brain, C6, 3T3, HeLa lysates using ENO3 antibody.



Immunocytochemistry analysis of ENO1/ENO2/ENO3(green) in HeLa using ENO1/ENO2/ENO3 antibody, and DAPI(blue)