


EPS15R Rabbit mAb

ABZ2940

Product Overview

Name	EPS15R Rabbit mAb
Accession(Primary)	Q9UBC2
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 94 kDa; Observed MW: 130 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:50

Target information (Q9UBC2)

Gene Name EPS15L1**About Gene Symbol**

Western blot analysis of EPS15R in K562 lysates using EPS15R antibody.

