

ERp72 Rabbit mAb


ABZ2942

Product Overview

Name	ERp72 Rabbit mAb
Accession(Primary)	P13667
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 73 kDa; Observed MW: 73 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:50-1:100

Target information (P13667)

Gene Name	PDIA4
About Gene Symbol	PDIA4 (Protein Disulfide Isomerase Family A Member 4) is a metabolic enzyme involved in cellular metabolism. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of ERp72 in PC-3 lysates
using ERp72 antibody.