

Ezrin Rabbit mAb


ABZ2944

Product Overview


Name	Ezrin Rabbit mAb
Accession(Primary)	P15311
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 69 kDa; Observed MW: 80 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:100-1:200; ICC/IF-1:200-1:500; FC-1:200-1:1000; IP-1:20-1:50

Target information (P15311)

Gene Name	EZR
About Gene Symbol	



Western blot analysis of Ezrin in K562, HeLa lysates using Ezrin antibody.



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using Ezrin antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.