

FES Rabbit mAb


ABZ2948

Product Overview

Name	FES Rabbit mAb
Accession(Primary)	P07332
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 93 kDa; Observed MW: 93 kDa
Recommended Dilution	WB-1:1000-1:5000; ICC/IF-1:100-1:200; FC-1:50-1:100

Target information (P07332)

Gene Name	FES
About Gene Symbol	



Western blot analysis of FES in Raji lysates using FES antibody.