


Fetuin A Rabbit mAb

ABZ2949

Product Overview

Name	Fetuin A Rabbit mAb
Accession(Primary)	P02765
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 39 kDa; Observed MW: 45-50 kDa
Recommended Dilution	WB-1:1000-1:5000

Target information (P02765)

Gene Name AHSG**About Gene Symbol**

Western blot analysis of Fetuin A in HepG2, 3T3, U937 lysates using Fetuin A antibody.

