


FGF1 Rabbit mAb

ABZ2950


Product Overview

Name	FGF1 Rabbit mAb
Accession(Primary)	P05230
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 17 kDa; Observed MW: 17 kDa
Recommended Dilution	WB-1:500-1:1000

Target information (P05230)

Gene Name FGF1**About Gene Symbol**

Western blot analysis of FGF1 in mouse heart lysates using FGF1 antibody.



Western blot analysis of FGF1 in C6 lysates using FGF1 antibody.