

FH/Fumarase Rabbit mAb


ABZ2952

Product Overview


Name	FH/Fumarase Rabbit mAb
Accession(Primary)	P07954
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 55 kDa; Observed MW: 49 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:200-1:500; ICC/IF-1:50-1:100; FC-1:200-1:500; IP-1:20-1:50

Target information (P07954)


Gene Name	FH
About Gene Symbol	FH (FumaRate HydRatase) is a metabolic enzyme involved in cellular metabolism. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.




Western blot analysis of FH/Fumarase in 3T3, Raw264.7 lysates using FH antibody.



Western blot analysis of FH/Fumarase in rat Brain, C6, CHO-K1, HeLa lysates using FH/Fumarase antibody.



Immunohistochemistry analysis of paraffin-embedded Human Cholangiocarcinoma using FH antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunocytochemistry analysis of FH/Fumarase (green) in Hela using FH/Fumarase antibody, and DAPI (blue).

