

Fibrillarin Rabbit mAb


ABZ2954

Product Overview

Name	Fibrillarin Rabbit mAb
Accession(Primary)	P22087
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 34 kDa; Observed MW: 34 kDa
Recommended Dilution	WB-1:1000-1:5000; ICC/IF-1:500-1:5000; FC-1:10-1:100

Target information (P22087)

Gene Name	FBL
About Gene Symbol	



Western blot analysis of Fibrillarin in K562, C6, 3T3 lysates using Fibrillarin antibody.