


Fibulin 1 Rabbit mAb

ABZ2955

Product Overview

Name	Fibulin 1 Rabbit mAb
Accession(Primary)	P23142
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 77 kDa; Observed MW: 77 kDa
Recommended Dilution	WB-1:1000-1:5000

Target information (P23142)

Gene Name FBLN1**About Gene Symbol**

Western blot analysis of Fibulin 1 in Human placenta, A549, 3T3 lysates using Fibulin 1 antibody.

