

Folate Binding Protein/FBP Rabbit mAb

ABZ2957


Product Overview

Name	Folate Binding Protein/FBP Rabbit mAb
Accession(Primary)	P15328
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 30 kDa; Observed MW: 40 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:200-1:1000


Target information (P15328)

Gene Name FOLR1

About Gene Symbol



Western blot analysis of Folate Binding Protein in HeLa lysates using Folate Binding Protein antibody.



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using Folate Binding Protein antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.