

Follistatin Rabbit mAb

ABZ2958


Product Overview

Name	Follistatin Rabbit mAb
Accession(Primary)	P19883
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Rat
Calculated MW	Calculated MW: 38 kDa; Observed MW: 38 kDa
Recommended Dilution	WB-1:500-1:1000

Target information (P19883)

Gene Name FST

About Gene Symbol



Western blot analysis of Follistatin in rat Brain lysates using Follistatin antibody.

