

FYB Rabbit mAb


ABZ2960

Product Overview

Name	FYB Rabbit mAb
Accession(Primary)	O15117
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 85 kDa; Observed MW: 130,120 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:200-1:500; FC-1:50-1:100

Target information (O15117)

Gene Name	FYB1
About Gene Symbol	



Western blot analysis of FYB in Jurkat lysates
using FYB antibody.