

Gelsolin Rabbit mAb


ABZ2965

Product Overview

Name	Gelsolin Rabbit mAb
Accession(Primary)	P06396
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 86 kDa; Observed MW: 86 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:100-1:200; IP-1:20-1:50

Target information (P06396)

Gene Name	GSN
About Gene Symbol	



Western blot analysis of Gelsolin in MCF-7, C6, HeLa lysates using Gelsolin antibody.