

GFPT2 Rabbit mAb


ABZ2967

Product Overview



Name	GFPT2 Rabbit mAb
Accession(Primary)	O94808
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 77 kDa; Observed MW: 77 kDa
Recommended Dilution	WB-1:500-1:1000; ICC/IF-1:200-1:500; FC-1:50-1:100; IP-1:20-1:50

Target information (O94808)


Gene Name	GFPT2
About Gene Symbol	



Immunocytochemistry analysis of GFPT2 (green) in HT-1080 using GFPT2 antibody, and DAPI (blue).



Western blot analysis of GFPT2 in rat spleen lysates using GFPT2 antibody.



Western blot analysis of GFPT2 in 3T3, HeLa
lysates using GFPT2 antibody