

Glucocorticoid Receptor Rabbit mAb

ABZ2969


Product Overview

Name	Glucocorticoid Receptor Rabbit mAb
Accession(Primary)	P06537
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 87 kDa; Observed MW: 94,91 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:200-1:2000; ICC/IF-1:200-1:500; FC-1:200-1:500


Target information (P06537)

Gene Name	Nr3C1
About Gene Symbol	


Rat skeletal muscle



Western blot analysis of Glucocorticoid Receptor in rat skeletal muscle lysates using Glucocorticoid Receptor antibody.



Western blot analysis of Glucocorticoid Receptor in K562, 3T3, HeLa lysates using Glucocorticoid Receptor antibody.



Immunocytochemistry analysis of Glucocorticoid Receptor(green) in Hela using Glucocorticoid Receptor antibody, and DAPI(blue)